


Shipping and Transport Instructions for ACL Units

To avoid any type of transport damage, we strongly recommend the use of ACL original packaging. If you do not have suitable transport packaging, we will be happy to send you a suitable one for your device - free of charge in case of a warranty claim, at extra charge outside the warranty. Improper packaging can lead to damage to individual components. For this reason, please ensure careful handling during transportation and loading. Always ensure gentle handling during charging and internal transport and avoid placing the unit on hard surfaces.

If you have any questions, please contact our service team at +49 341 23078-60 or by email at service@acl.de.

OR-PC®/ OR-MD® 15

1. Remove all cables, mounting systems and other accessories.
2. Place the foam in front of you and bend the following foam parts into the respective opening.



3. Turn the foam around, pull it slightly apart and place the device with the backside into the foam.
4. Take the corners of the two foam parts and clamp it over the display on the opposite side.



5. Hold the box placed upright and slide the device including the foam parts in the provided transport box.





6. Store the accessories (power cable, etc.) in the accessory box. Place the accessory box on the upper part as shown on the picture and close the box.



Shipping and Transport Instructions for ACL Units

To avoid any type of transport damage, we strongly recommend the use of ACL original packaging. If you do not have suitable transport packaging, we will be happy to send you a suitable one for your device - free of charge in case of a warranty claim, at extra charge outside the warranty. Improper packaging can lead to damage to individual components. For this reason, please ensure careful handling during transportation and loading. Always ensure gentle handling during charging and internal transport and avoid placing the unit on hard surfaces.

If you have any questions, please contact our service team at +49 341 23078-60 or by email at service@acl.de.

OR-PC®, OR-MD® 21	OR-PC®, OR-MD® 24
<ol style="list-style-type: none"> 1. Remove all cables, mounting systems and other accessories. 2. Place the unit between the black foam parts. 	<ol style="list-style-type: none"> 1. Remove all cables, mounting systems and other accessories. 2. Place the unit between the black foam parts. Remove the edge and throw away the remaining parts. 

3. Hold the box placed upright and slide the device including the foam parts in the provided transport box.





4. Store the accessories (power cable, etc.) in the accessory box. Place the accessory box on the upper part as shown on the picture and close the box.



Shipping and Transport Instructions for ACL Units

To avoid any type of transport damage, we strongly recommend the use of ACL original packaging. If you do not have suitable transport packaging, we will be happy to send you a suitable one for your device - free of charge in case of a warranty claim, at extra charge outside the warranty. Improper packaging can lead to damage to individual components. For this reason, please ensure careful handling during transportation and loading. Always ensure gentle handling during charging and internal transport and avoid placing the unit on hard surfaces.

If you have any questions, please contact our service team at +49 341 23078-60 or by email at service@acl.de.

OR-PC®/ OR-MD® 19	OR-PC®/ OR-MD® 27, 32, 42, 43, 55
<ol style="list-style-type: none"> 1. Remove all cables, mounting systems and other accessories. 2. Place the unit between the plastic foils. Press both foils onto each other. The unit shall not have any space to move between the foils as shown on the picture below:  3. Place the accessory box on the upper plastic foil. Use a knife to cut down all four carton corners until the foils have no room left for movements inside the box as shown on the picture below:  <p>Close the box after the proper placement of the unit according to the instructions above.</p> 	<ol style="list-style-type: none"> 1. Remove all cables, mounting systems and other accessories. 2. Attach the airbags to two opposing sides of the unit and place it in the carton. <p>The carton shall be fixed with tension straps onto a pallet in an upright position. Transporting the carton in a lying position might damage the unit!</p> <p>The box must not be larger than the pallet as shown on the picture below:</p> 